



CryoCell

CRYOCELL INSERT FOR THE ESL PLATFORM

ULTRA LOW ABLATION TEMPERATURES
FROM THE SIMPLEST DESIGN

FEATURES

- CLEAN PELTIER DESIGN
- EFFICIENT TCUBE BENCHTOP CHILLER
- TWOVOL2 AND TWOVOL3 CAPABILITY
- ACTIVEVIEW2 SOFTWARE ADD-IN

High performance cooling capabilities
with 1 °C set point intervals

Temperature stability of 1 °C over extended
periods of time

Quick cooling of sample space to set point

Efficient switch to and from CryoCell applications



© Elemental Scientific Lasers LLC | 685 Old Buffalo Trail | Bozeman, MT, 59715 | United States
Tel: + 406 586 3159 | lasers@icpms.com | www.icpmslasers.com