

ESL266macro

**LOWEST DETECTION LIMITS FOR
BULK ANALYSIS APPLICATIONS**

ACHIEVE LOWEST DETECTION LIMITS FOR YOUR
BULK ANALYSIS APPLICATION WITH ICP AND ICPMS

FEATURES

- LARGE SPOT SIZE
- AUTOMATION
- ACTIVEVIEW2 SOFTWARE


Up to 1 mm spot size for improved sensitivity of bulk and spatial analysis

100 mm x 100 mm sample chamber for large samples and hands-off automation – compatible with robotic sample handling

Intuitive, feature-rich laser ablation control software – designed around your workflows

Multi-layered viewing for outstanding sample navigation and quick experimental setup

Improved ICP/ICPMS compatibility



© Elemental Scientific Lasers LLC | 685 Old Buffalo Trail | Bozeman, MT, 59715 | United States
Tel: + 406 586 3159 | lasers@icpms.com | www.icpmslasers.com