


TwoVol2

UNSURPASSED POSITIONAL REPRODUCIBILITY AND STAGE ACCURACY

THE HIGHEST PERFORMING LASER ABLATION CHAMBER ON THE MARKET



FEATURES

- TYPHOON PURGE
- STAGE RETURN ACCURACY $<1.5 \mu\text{M}$
- COMPATIBLE WITH ESL193, ESL213, ESLFEMTO, AND ESLIMAGE

Efficient Typhoon purge mechanism enables stable operation after just 5 minutes and spatial reproducibility of $<2\%$ RSD

Non-cantilevered stages to accurately target features with minimal wear

Leak free operation with no complex corner seal and accessible tubing for easy replacement



© Elemental Scientific Lasers LLC | 685 Old Buffalo Trail | Bozeman, MT, 59715 | United States
Tel: + 406 586 3159 | lasers@icpms.com | www.icpmslasers.com